

exchange[®]

introducing

mysafe



ensuring safe disposal of 1ml
insulin type syringes

www.harmreduction.co.uk



The sharing of used injecting equipment makes people who inject drugs particularly vulnerable to HIV and Hepatitis C.¹

1ml insulin type syringes require a method of ensuring safe disposal after use. MySafe is designed to lock onto the outer barrel of the syringe, securing the needle within the MySafe. The cannula is bent upon entry rendering it unusable.

how it works



Step 1

Ready the MySafe for use after injecting



Step 2

Insert used syringe into the MySafe



Step 3

Push until the syringe locks into the MySafe

get in touch

Visit: www.harmreduction.co.uk

Email: info@harmreduction.co.uk

Call: **01495 235 800**

Twitter: [@FMG_Exchange](https://twitter.com/FMG_Exchange)

¹ World Drug Report 2014, UNODC (Available at <http://www.unodc.org/wdr2014>)

